

FUMEHOOD SURVEY

HOOD ID# _____

FACE VELOCITY _____ **FPM**

NAME	DATE	RESULT

GOOD WORK PRACTICES FOR FUMEHOOD USE

- Store chemicals and equipment outside of hood
- Keep containers and equipment >6" from front/back of hood
- Provide space under large equipment for proper air flow
- Keep sash closed or at lowest possible usable position
- Do not conduct perchloric acid digestion in this hood
- Prevent release of chemicals into fumehood sink
- Open chemical storage containers only when dispensing
- Avoid open container procedures when using volatiles

